

Balance SL-16

Liquid 16oz

- **Broad-spectrum antimicrobial**
- **Anti-inflammatory (as inflammation is often caused by bacteria, virus, or mold)**
- **Non-toxic, even if used multiple times per day for weeks Does not encourage resistant microbial strains**
- **Immune system support**
- **Does not interfere with antibiotics**

Balance SL-16 is a highly effective antimicrobial preparation composed of pure silver complexed with purified water. Balance SL-16 uses Silver Sol Technology®, the term “sol” being a designation of a mineral permanently distributed into the structure of water. The silver becomes a part of that water molecule permanently, so it will not fall out of solution or suspension.

Balance SL-16 contains 15 ppm of purified silver per serving. This is an incredibly powerful, non-toxic form of silver, with zero build-up in the body, so it does not cause argyria (blue man syndrome).

Made with non-GMO ingredients.

Balance SL-16 is a highly effective antimicrobial composed of pure silver complexed with purified water. Balance SL-16 utilizes a technological breakthrough allowing for the production of a silver particle small enough to enter a single red blood cell,

enabling the silver to travel through capillaries, the tiniest blood vessels in the body, and inactivate viruses at the cellular level. Also, Balance SL-16 uses technology that permanently distributes the mineral into the structure of the water, so it will not fall out of solution or suspension. As a result, the pure silver in Balance SL-16 has been shown to provide a whole host of health benefits, from killing drug-resistant bacterial strains to inhibiting viral replication as well as stimulating stem cell production. It can be safely used both orally and topically with no known side effects.

Balance SL-16 – The superior silver; distinct from traditional colloidal (ionic) silvers Balance SL-16 is completely distinct from traditional colloidal (ionic) silvers of the past, with many important chemical, physical and vibrational differences separating it far beyond other forms available. In fact, while no other liquid silver antimicrobial product had received a patent for over 80 years, a US patent was issued in 2006 for this silver composition.

Colloidal silver works by chemical action, requiring direct contact with microbes, where each colloidal silver molecule can only kill one bacteria. Also, colloidal silver is simply a suspension of a silver molecule and an oxygen molecule combined together. Particles in a suspension, while too small to be seen, are large enough to settle to the bottom over time. Here, the silver components can fall out of solution,

causing them to aggregate in such areas as the blood, fat, and the brain. Not surprisingly, silver suspensions contain up to 300,000 ppm of silver compared to Balance SL-16 which is effective at a much lower dosage of only 5–30 ppm. Argyria (blue man syndrome, where silver aggregates in the body) has only been seen in cases of colloidal silver use after extended high dosing.

Three main mechanisms of action for Balance SL-16

1. Catalytic Action – far more effective than chemical action (how colloidal silver works) Balance SL-16 does much more than chemical action because it works by what is known as catalytic action. Balance SL-16 is made up of four atoms of oxygen and four atoms of silver (Ag_4O_4). Here, one silver atom is missing three electrons while another is only missing one. The molecule tries to balance itself out, so a single, high-energy particle goes back and forth

at high speed between the two atoms, which consequently gives it catalytic capabilities. Hence, one electron destroys the bacteria (or other pathogen), while the other electron is being recharged. Thus, this catalytic action allows the Balance SL-16 molecule to destroy over and over, again and again. The result is that Balance SL-16 can eradicate thousands of times more pathogens than a colloidal or ionic silver, allowing it to be much more effective at much lower, safer concentrations.

2. Resonance This high energy that is constantly occurring creates a resonant frequency, where the silver is resonating at a frequency that is selectively destructive to pathogens. This energy has been measured at 890-910 terahertz, which is the same frequency at which germicidal ultraviolet light resonates. Again, the silver particles in Balance SL-16 are tiny enough to be absorbed by red blood cells and they resonate at the perfect frequency to destroy pathogens.

3. Viral Disruption Balance SL-16 can bind with the charged, genetically incomplete viral molecules, preventing viral replication. And because healthy human DNA does not contain a magnetic charge, it remains unaffected by the charge of the Balance SL-16.

THIS INFORMATION IS PROVIDED FOR THE USE OF PHYSICIANS AND OTHER LICENSED HEALTH CARE PRACTITIONERS ONLY. THIS INFORMATION IS INTENDED FOR PHYSICIANS AND OTHER LICENSED HEALTH CARE PROVIDERS TO USE AS A BASIS FOR DETERMINING WHETHER OR NOT TO RECOMMEND THESE PRODUCTS TO THEIR PATIENTS. THIS MEDICAL AND SCIENTIFIC INFORMATION IS NOT FOR USE BY CONSUMERS. THE DIETARY SUPPLEMENT PRODUCTS OFFERED BY DESIGNS FOR HEALTH ARE NOT INTENDED FOR USE BY CONSUMERS AS A MEANS TO CURE, TREAT, PREVENT, DIAGNOSE, OR MITIGATE ANY DISEASE OR OTHER MEDICAL CONDITION.

Safety

After using Balance SL-16, 99% of the silver leaves the body in 24 hours; all is cleared by 48 hours. This occurs because the silver in Balance SL-16 has been structured into the water, allowing the water molecule to get absorbed into the blood stream and pass through circulation. It ultimately gets excreted unchanged, which means no harmful metabolites are formed.

A study on humans with the AIDS virus demonstrates Silvercillin's safety and beneficial properties in those who are immune-compromised. Dose for four months. Results showed a normalization of T-lymphocyte counts as well as restoration of normal or close-to-normal body weights.

Does not harm good bacteria

Balance SL-16 is water-soluble, rendering it unable to penetrate fats. Because of this, it cannot harm good bacteria, such as lactobacillus, which include a milk-like fat surrounding the bacteria. This difference sets them apart from most other bacteria molecules. Healthy cells in the body also contain a protective lipid membrane (unlike the thin cell walls of pathogens) with a balanced charge, making them well protected from the silver's oxide.

Safe with Antibiotics. Balance SL-16 is not only safe to take with antibiotics, but can actually work synergistically to increase the effectiveness of specific antibiotics. With its different mechanism of action than that of antibiotics, Balance SL-16 will help to “mop up” the resistant pathogens that the antibiotic may have left behind. A 2006 study published in Current Science looked at 19 common antibiotics in combination with a silver solution against seven different microbial organisms. The results of this study showed the increased efficacy of some antibiotics up to 10-fold, as the antibiotic-silver combination allowed for a more complete clearing of the pathological organism.

Anti-Bacterial

The active ingredient in Balance SL-16 has been shown to inhibit and/or kill an extensive array of harmful bacteria including MRSA (methicillin resistant Staphylococcus aureus), which can be deadly. A Brigham Young University report found this silver to be effective at killing every pathogenic bacterium it was tested against at less than 10 ppm. Among this list were salmonella and E. coli (food poisoning), Klebsiella pneumonia and oxytoca (pneumonia), Haemophilus influenza (ear infections, meningitis), Enterobacter aerogenes (UTIs), and Streptococcus mutans and gordonii (tooth decay).

Anti-Viral (via magnetic disruption)

Balance SL-16 technology has the ability to inactivate viruses at the DNA/RNA level. Incomplete segments of DNA/RNA are found in viruses and carry a slight magnetic charge. The silver in Balance SL-16 acts like a magnet, attracting these charged particles and binding the DNA so tightly that it inactivates the virus and prevents virus replication. This same mechanism of action does not take place with healthy cells, as their thicker, more protected cell membranes carry a balanced charge, which protects them from the silver’s magnetic attraction.

Anti-Yeast, Mold, Fungi

Fungal infections are often refractory which can make them drug-resistant and difficult to treat. Numerous studies show that Silvercillin’s lethal properties are

also evident in regard to mold, fungi, and yeast, making it an excellent choice for use against Candida albicans. A double-blind study revealed that silver was effective at inhibiting the growth of and killing Candida albicans at a dose of 10 ppm after only two hours of exposure. Also consider using Balance SL-16 in the oral cavity to help prevent fungal infections (Candida/thrush) and oral plaque build-up.

Manufactured for Food Intelligence 22311 Stevens Creek Blvd Cupertino, CA 95014 www.mor-nutrition4life.com

Recommended Use: As a dietary supplement, take 5 mL (approx. one teaspoon) orally per day, or as directed by your health care practitioner.

Supplement Facts

Serving Size 5 mL (approx. 1 teaspoon)

Servings Per Container about 95

Amount Per Serving	% Daily Value
Purified Silver	75 mcg *

*Daily Value not established.

Other Ingredients: Purified water.